

Saving Results

Cytoscape provides a number of ways to save results and visualizations:

- As a session: **File → Save, File → Save As...**
- As an image: **File → Export → Network to Image...**
- To the web: **File → Export → Network to Web Page...** ([Example](#))
- To a public repository: **File → Export → Network to NDEx**, or **File → Export → Collection to NDEx**
- As a graph format file: **File → Export → Network to File.**

Formats:

- CX JSON
- Cytoscape.js JSON
- GraphML
- PSI-MI
- XGMML
- SIF

